

PANBIO™ COVID-19 Ag RAPID TEST DEVICE

SCALE UP COVID-19 ANTIGEN TESTING

For patients suspected of current COVID-19 infection.

Requires no instrumentation and provides results in 15 minutes, making it a valuable tool for mass testing in decentralized settings.

Accessible, portable and scalable option for COVID-19 testing.

May also be useful for supporting public health strategies, such as contact tracing and large-scale testing of people suspected of having an active infection.

Product not available in all countries. Not approved for sale in the United States.



CONTACT SALES

TECHNICAL SUPPORT

CHAPTER NAVIGATION



BENEFITS

- Rapid answers for patients: test results in 15 minutes
 - Reliable, high performance*: Sensitivity: 91.4% (94.1% for samples with Ct values ≤ 33) / Specificity: 99.8%
 - Deploy at point of care in a variety of settings
 - Accessible solution enables large scale testing
 - Self-contained tube with "break off" swab minimizes staff exposure
- We use cookies, just to track visits to our website, we store no personal details.

ACCEPT

[PRIVACY POLICY](#)

SPECIFICATIONS



- ✓ Test time: 15 minutes
- ✓ Sensitivity*: 91.4% (94.1% for samples with Ct values ≤ 33)
- ✓ Specificity*: 99.8%
- ✓ Storage: 2°C–30°C
- ✓ CE MARK
- ✓ Sample type: Nasopharyngeal swab
- ✓ Quantity: 25 tests/kit

PRODUCT CODE



VIDEOS

We use cookies, just to track visits to our website, we store no personal details.

ACCEPT

[PRIVACY POLICY](#)

Panbio™ COVID-19 Ag Rapid Test Device

NASOPHARYNGEAL SWAB TEST PROCEDURE



PANBIO™ COVID-19 AG RAPID TEST DEVICE NASOPHARYNGEAL SWAB TEST PROCEDURE

RELATED PRODUCTS

Learn how to use the Panbio™ COVID-19 Ag Rapid Test.



*Cassette design may vary

PANBIO™ COVID-19 IgG/IgM RAPID TEST DEVICE

We use cookies, just to track visits to our website, we store no personal details.

ACCEPT

[PRIVACY POLICY](#)

REFERENCES



+

STAY INFORMED

Sign up to receive valuable updates from Abbott.

THE SYMPHEOS™ SYSTEM

[CLICK HERE TO SIGNUP](#)

NEED TECHNICAL SUPPORT?

Click here to get in touch with our technical support team.

TECHNICAL SUPPORT

We use cookies, just to track visits to our website, we store no personal details.

ACCEPT

PRIVACY POLICY